



# Accessories

## Adjustable monitor stand

Specially designed to support bedhead monitors, allowing its orientation on a vertical axis and its inclination by means of a lower control, it enables visualization from any point.

- The rail version is attached by means of a clamp with two knobs, while the wall version is attached thanks to a slide vertical profile. The Ø38 mm version is fixed by two flanges.
- Every model has a practical basket coated in white.
- The standard version includes a stand with buckles. In addition there are specific versions adjusted to *General Electric*, *PHILIPS SureSigns* and *MP* monitors, as well as some other brands on the market.



### TECHNICAL SPECIFICATIONS

Manufactured in anodized aluminium (profiles) and zinc steel.

- The standard platform is 300 mm wide by 180 mm deep. It includes a buckle to secure the monitor.
- The horizontal arm is 220 mm long and has an inferior reinforcement.
- The vertical profile is 500 mm long and allows an adjustment of 270 mm within the total height of the set, thanks to its rails.
- **Rotations:**
  - Lateral rotation: 150°.
  - Tilt from 0° to 15°.
- **Weight:**
  - Rail version: 2,7 kgr.
  - Wall version: 3,3 kgr.
  - Ø38mm tube version:
- **Maximum load:** 25 kg.