



# Accessories

## Ambulance equipment



This equipment is fit for oxygen supply and suction. It can contain an auxiliary oxygen outlet nozzle, which is fed from the same source.

- It is especially designed to be used in ambulances, but it can be adapted to any other place with an 4.5 bar oxygen source.
- There are multiple configurations, in order to adjust to concrete needs.
- Available in wall and rail version with different equipment configurations.

### TECHNICAL SPECIFICATIONS

- Dial flowmeter or ball flowmeter, Venturi-effect suction and an additional quick outlet can be configured by the client..
- Supply pressure: 4.5 bar.