

CLAIRE

Decorative horizontal bedhead

The **CLAIRE bedhead unit** is a hospital bedhead designed to respond to the highest aesthetics and functionality requirements, becoming another architectural element of the environment. The mechanisms, except the gas outlets, are incorporated to its lower profile, parallel to the ground, thus being hidden by the front cover.

Starting on some basic components, our technical office will adapt the dimensions, colours and final equipment to the specific needs of each project upon demand.



TECHNICAL SPECIFICATIONS

Manufactured in:

- Structure and covers manufactured in extruded aluminium, forming three independent and removable compartments.
- Translucent covers in polycarbonate.
- Covers ends in ABS.

Its dimensions are 122 mm depth and 233 mm height.

Finishings:

- Main body lacquered: RAL 9010.
- Covers lacquered in the following colours:
 - Red: RAL 3003
 - Blue: RAL 5015
 - Green: RAL 6019
 - Grey: RAL 7035
 - White: RAL 9010
- The front cover can be customized with different bedhead finishings for different services areas, and the image of the hospital can be incorporated.

In compliance with UNE-EN ISO 11197 standard.