



Accessories

Collection jars

The Collection Jar is designed to collect suction fluids. It incorporates an overflow trap and is available in different volumes: 300, 500, 1000 and 2000 ml.

- Reusable after sterilization thanks to its metallic connection and the materials used in the manufacture process.



TECHNICAL SPECIFICATIONS

300 cc Jar:

- Manufactured in *PC, PA, chrome-plated brass and NBR.*
- Connection to the vacuum regulator with a DISS 3/4-16 UNF nozzle.
- Sterilizable in autoclave up to 121 °C and 1,3 kPa (max. 40 cycles)
- Maximum level:
 - 11VSR0300: 300 cc.



500 cc, 1000 cc and 2000 cc Jars:

- Manufactured in *PC, PA y lacquered steel (RAL 9010).*
- Connection to the vacuum regulator by flexible hose.
- Sterilizable in autoclave up to 120 °C.
- Maximum level:
 - 11VSR0500: 500 cc.
 - 11VSR1000: 1000 cc.
 - 11VSR2000: 2000 cc.

