

Outlet converter



JUST PLUG C is especially useful when the equipment connectors are different in standards from those of the valves..

- The set is composed by a connector and an outlet for the same gas but in compliance with different standards.

TECHNICAL SPECIFICATIONS

- Manufactured with Chromed-plated brass, stainless steel and NBR.
- In compliance with:
 - NF S 90-116 standard.
 - DIN 13260-2 standard.
 - BS 5682 standard.
 - SS 87 524 30 standard.
 - CGA V-5 (DISS) standard.