




# Accessories

## Drawers

Set of drawers attachable to a rail by means of two pushbutton clamps.

- Composed by an upper tray and drawers in different formats.
- The inner buckets can be easily be extracted out of the drawer.
- Thanks to its slide partitions, it is possible to establish different divisions inside the drawers.



### TECHNICAL SPECIFICATIONS

Infusion pump holders for rail (*ref.: 13SPBM0000*):

- Manufactured in anodized aluminium of Ø25x1180 mm dimensions with a telescopic extension for four bags of serum.
- Anchored to two rails by means of pushbutton clamps.

Infusion pumps holders for rail (*ref.: 13SCCG3500*):

- Manufactured in anodized aluminium of Ø25x350 mm dimensions with a telescopic extension for two bags of serum.
- Anchored to a rail by means of a clamp with knob.