

Outlet duplicators



JUST PLUG D is especially designed to achieve two outlets for the same gas, in systems where there is only one.

- The set is composed by a connector, a distributor and two outlets for the same gas in compliance with the same or a different standard.

TECHNICAL SPECIFICATIONS

- Manufactured with Chromed-plated brass, stainless steel and NBR.
- In compliance with:
 - NF S 90-116 standard.
 - DIN 13260-2 standard.
 - BS 5682 standard.
 - SS 87 524 30 standard.
 - CGA V-5 (DISS) standard.