

# JUST PLUG Series

European Standards Connectors and Outlets

Ref.:

**JUST PLUG P** connectors with thread for equipment or for a flexible hose.

- Available for Oxygen, Vacuum, Medical Air and Nitrous Oxide (*please consult our staff about other gases*).
- BSP 1/4" x 19 thread for vacuum and BSP 1/8" x 28 for the other gases.



## TECHNICAL SPECIFICATIONS

Manufactured standards:

- Standard **NF S90-116**.
- Standard **DIN 13260**.
- Standard **BS 5682**.
- Standard **SS 87 524 30**.
- Standard **NIST BS 5682**.
- Standard **UNE-EN ISO 9170-2 Type 1L SEGA**.

Manufactured in *chrome-plated brass*.

**C/ Tungsteno 11-15.**  
**Arganda del Rey, 28500, Madrid**  
**Tel.:+34 91 871 99 50**  
[www.itcsal.com](http://www.itcsal.com)      [itcsal@itcsal.com](mailto:itcsal@itcsal.com)

09/2017-Rev.02



# JUST PLUG Series

European Standards Connectors and Outlets

Ref.:

**JUST PLUG T** outlets for bedheads and ambulances, either with a wall box or embedded to the wall.

- Available for Oxygen, Vacuum, Medical Air and Nitrous Oxide (*please consult our staff about other gases*).
- With  $\varnothing 10\text{mm}$  copper tube to weld the outlets.



## TECHNICAL SPECIFICATIONS

Manufactured standards:

- Standard **NF S90-116**.
- Standard **DIN 13260**.
- Standard **BS 5682**.
- Standard **SS 87 524 30**.
- Standard **NIST BS 5682**.
- Standard **UNE-EN ISO 9170-2 Type 1L SEGA**.

Manufactured in *chrome-plated brass*.

**C/ Tungsteno 11-15.**  
**Arganda del Rey, 28500, Madrid**  
**Telf.: +34 91 871 99 50**  
[www.itcsal.com](http://www.itcsal.com)      [itcsal@itcsal.com](mailto:itcsal@itcsal.com)

09/2017-Rev.02

