



Flexible hoses designed for the conduction of medical gases from the pipe network outlets to the points of consumption, such as anaesthesia devices, ventilators, etc...

- Available for Oxygen (white), Vacuum (yellow), Medical Air (black and white), Nitrous Oxide (blue), Carbon dioxide (grey) and Nitrogen (black).

TECHNICAL SPECIFICATIONS

- Composed by a PVC interior pipe, a middle polyester fibre reinforcement and an exterior PVC cover.
- Antistatic hose.
- **Maximum test pressure:** 15 bar.
- **Maximum work pressure:** 10 bar.
- **Dimensions:**
 - Vacuum: inner Ø8 mm and outer Ø14 mm.
 - Other gases: inner Ø6 mm and outer Ø12 mm.
- Branded with: **LOT, GAS and ISO 5359**
- In compliance with **UNE-EN ISO 5359** standard