



Accessories

Flow selector



The **Flow selector** is a device able to select, thanks to different flowmeters, a dry or a humid outlet (the latter through a humidifier bottle), without dismantling its components:

- The **Flow selector** is inserted between the flowmeter and the humidifier bottle. For this, they must have a standard DISS connection.
- It has a lower outlet for the DISS humidifier bottle (not included) and a side direct outlet, easily selected by means of a quarter of a turn knob.

TECHNICAL SPECIFICATIONS

- Manufactured in natural *anodized aluminium, chrome-plated brass and nylon*.
- DISS 9/16" connection.
- **Dimensions:** Ø25 x 60 mm.
- **Weight:** 100 g.