



# Accessories

Thermic blanket ~ Guedell Airways ~ Jaw spreader ~ Extension lead ~ Tongue clamp



Emergency

- **Thermic blanket** for protection against cold, heat and humidity. Presented individually in a transparent plastic bag with instructions for use.
- **Guedell airways**, airways tubes easily identifiable through colour codes.
- **Helicoidal jaw spreader** to allow the opening of air ways.
- 2m long **extension lead**, which facilitates an oxygen flow from a flowmeter to a manual resuscitator.
- **Tongue clamp**, to avoid patients from choking on their tongues when unconscious.

## TECHNICAL SPECIFICATIONS

### Thermic blanket:

- The kit is composed of: 160 x 210 cm sheet manufactured in polyester metalized with aluminium on one side and lacquered golden on the other.

### Guedell airways:

- 09TUGU1000 Tube n° 1 of 70 mm with white connector.
- 09TUGU2000 Tube n° 2 of 80 mm with green connector.
- 09TUGU3000 Tube n° 3 of 110 mm with orange connector.

### Jaw spreader:

- Manufactured in Polycarbonate.

### Extension lead:

- Manufactured in PVC.

### Tongue clamp:

- Manufactured in Polyamide.