



# Accessories

## NEXUS cart - Metallic Shelf

Shelves attached to different cart and bedhead structures.

- Designed with a small projection along the perimeter to prevent objects from falling.
- They have two side flanges in order to adapt to the cart structures, enabling a quick regulation.



### TECHNICAL SPECIFICATIONS

#### Manufactured in:

- Shelf manufactured with lacquered steel sheet (RAL 9010).
- The finishing elements are made of ABS in different RAL colours.

#### Finishing:

- Lacquered: **RAL 9010**.
- Finishing elements available in the following colours:
  - Red: **RAL 3003**
  - Blue: **RAL 5015**
  - Green: **RAL 6019**
  - Grey: **RAL 7035**

Maximum load: 25 kg.