



NEXUS LT

Ceiling-pendant

The ceiling-pendant NEXUS LT are designed to integrate the supply of medical gases and electricity and communications plugs, and are specially indicated for surgical areas and intensive care and neonatal units.

- The ceiling-pendant comprises a support to the ceiling and a horizontal profile containing six isolated compartments (three in the LT-S version), which can be equipped differently depending on the needs of each project.
- Its design displays the gas outlets in working arrangement, allowing the operator to connect flowmeters and vacuum regulators directly during delicate situations, avoiding connection by flexible pipes.
- Complemented with rotating racks of different heights, trays, drawers and diverse accessories from our rail commercial line.



Bedhead units

TECHNICAL SPECIFICATIONS

Manufactured in:

- Structure and covers in extruded aluminium, forming independent removable compartments.
- Translucent covers in polycarbonate.
- Covers ends in polystyrene.

Its dimensions are 414 mm depth by 255 mm height and 500 mm length.

Finishings:

- Main body lacquered: **RAL 9010**.
- Covers lacquered in the following colours:
 - Red: **RAL 3003**
 - Blue: **RAL 5015**
 - Green: **RAL 6019**
 - Grey: **RAL 7035**
 - White: **RAL 9010**

In compliance with UNE-EN ISO 11197 standard.