



# Accessories

## PUSH System accessories

The PUSH System of the JUST VAC suction units allows a quick and easy replacement of its accessories. This is specially designed to protect the device from possible overflowing bottles or from particle accumulation towards its interior.

- Pressing the **PUSH** button, the accessory is released for its replacement.
- In case the accessories do not incorporate the PUSH System, we also have devices with 1/4" and 1/8" BSP nozzles.



### TECHNICAL SPECIFICATIONS

- Disposable filters made from *Polystyrene and Polyurethane Foam*.
- Connectors in *chrome-plated brass*.
- Devices PUSH in *ABS, chrome-plated brass and NBR*.