



Accessories

Recognition LED lamp

The newest recognition lamp is the lamp of cold light, which relies on LED diodes technology. It is indicated for any medical specialty, but mainly for those places where a total absence of calorific radiation is required and a very good light quality is still needed. Note the long life of the LED diodes, almost unlimited, and its low energetic consumption, lower than 10 W.

Its design and robustness make it fit for every kind of use in clinics, hospitals and consultation facilities. Rail version and Ø38mm tube version available, thus being able to adapt to any work area. High mobility thanks to the high quality and reliability of the flex tube. Together with the lozenge next to the cupola, they make this device able to adapt to an infinity of positions.



TECHNICAL SPECIFICATIONS

- Voltage: 95-240V, 50-60 Hz.
- Light source: LED.
- Luminous angle: 6°.
- Power: 10 W
- Short circuit protection: Yes.
- Thermal protection: Yes.
- Lifetime: 50.000 hours.
- Luminous intensity: 42 Klux / 0,5 m; 12Klux / 1,0 m
- Colour temperature: 4.750 K
- Light field diameter to 1m: 12 cm