



# Accessories

## NEXUS twistable cart / twistable and sliding cart

This cart is capable of turning gently on a vertical axis and moving on the horizontal axis, holding equipment and accessories on the patient's bedhead, without the need of supports on the floor or on the walls.

- There are two standard lengths adaptable to each need. Both racks have four lateral profiles and two inferior ones for the anchoring of various accessories from our rail commercial line.
- Accessories such as shelves, drawers and auxiliary rails can be easily accommodated to vertical tubes with removable flanges, allowing a quick change in the configuration of the cart.



Bedhead units

### TECHNICAL SPECIFICATIONS

#### Manufactured in:

- Manufactured with lacquered steel profiles and  $\varnothing 38$  mm stainless steel tubes.
- They have a frontal and an inferior system that enables displacement and turn blocking.
- Available in two widths and two standard heights, although the latter can be adapted to the client's needs.

#### Finishing:

- Lacquered: **RAL 9010**.
- Finishing elements available in the following colours:
  - Red: **RAL 3003**
  - Blue: **RAL 5015**
  - Green: **RAL 6019**
  - Grey: **RAL 7035**

Maximum load: 150 kg.