



Accessories

Shelves

Shelves attached to the rail by two pushbutton clamps.

- Designed with a little projection along the perimeter to prevent objects from falling.
- They have two regulators in the wall holders that guarantee a perfect horizontal levelling.



TECHNICAL SPECIFICATIONS

Manufactured:

- Tray manufactured in lacquered steel (RAL 9010).
- Finishing elements are made of ABS in different RAL colours.

Finishing:

- Lacquered: **RAL 9010**.
- Finishing elements available in:
 - Red: **RAL 3003**
 - Blue: **RAL 5015**
 - Green: **RAL 6019**
 - Grey: **RAL 7035**

Maximum load: 15 kgr.