



Stretchers and stands

Stretcher stands



Bellow' stand. Tool-free opening and closing. Blockage performed by the own stretcher.


TECHNICAL SPECIFICATIONS

Manufactured with anodized aluminium frames.

Dimensions:

- Closed: 60 x 19 x 88 cm
- Open: 60 x 114 x 88 cm
- Weight: 6,8 Kg

Pre-assembled stand. Set composed by two half stands that must be assembled before their.



TECHNICAL SPECIFICATIONS

Manufactured with anodized aluminium frames.

Dimensions:

- Closed: 6 x 59 x 144 cm
- Open: 120 x 59 x 78 cm
- Weight: 6,25 Kg



Foldable stand. A pin allows the stretcher to be folded until all components are placed on the same level.

TECHNICAL SPECIFICATIONS

Manufactured with anodized aluminium frame.

Dimensions:

- Closed: 6 x 57 x 183 cm
- Open: 102 x 57 x 64 cm
- Weight: 9,0 Kg