



Accessories

Steerable rack for rail

Columns designed for patients with special needs, to be placed at the bedhead.

- Attached to both an upper and a lower rail by means of four clamps with knobs.
- It is possible to configure them either with three trays, or with two trays and a drawer, being these in different formats.
- The whole set can be turned horizontally and both the trays and the drawers can be moved vertically.



TECHNICAL SPECIFICATIONS

Manufactured in:

- Structures manufactured in anodized aluminium, with Ø38 mm polish stainless steel tubes.
- Tray manufactured in lacquered steel.
- Drawers manufactured in aluminium with finishing elements in ABS.

Finishing:

- Lacquered in white: **RAL 9010**.
- Finishing elements also available in:
 - Red: **RAL 3003**
 - Blue: **RAL 5015**
 - Green: **RAL 6019**
 - Grey: **RAL 7035**

Maximum load: 60 kgr.