



Accessories

Supports for droppers

Hangers of different types for serum bags to prevent obstacles on the floor.

- The mount incorporates a clamp with a knob, which anchors the hanger firmly to the rail. In addition, a knob allows vertical repositioning for an easier access to the hooks.
- There is a $\varnothing 38,0$ mm bar mounting option using flanges.



TECHNICAL SPECIFICATIONS

- Manufactured in *anodized aluminium*.
- Bars of either $\varnothing 13 \times 450 \times 850$ mm (with two folds) or $\varnothing 13 \times 450 \times 1000$ mm (in right angle) dimensions.