



Accessories

Infusion pumps holder

Supports specially designed for the anchoring of perfusion pumps in the patient's bedside without the need to rest on the floor, both in its rail version and its Ø38.0 mm bar version.



TECHNICAL SPECIFICATIONS

Infusion pump holders for rail (*ref.: 13SPBM0000*):

- Manufactured in anodized aluminium of Ø25x1180 mm dimensions with a telescopic extension for four bags of serum.
- Anchored to two rails by means of pushbutton clamps.

Infusion pumps holders for rail (*ref.: 13SCCG3500*):

- Manufactured in anodized aluminium of Ø25x350 mm dimensions with a telescopic extension for two bags of serum.
- Anchored to a rail by means of a clamp with knob.