



TECHNICAL WAY

Basic horizontal bedhead

The **TECHNICAL WAY Bedhead unit** is a wall bedhead that integrates the supply of medical gases, power, telephone and computer plugs, etc...

It can be mounted either in single, double or triple bars, effectively separating the electricity cables and the gas pipes. Our technical office will adjust the dimensions, colours and final equipment to the specific needs of each project upon demand.



TECHNICAL SPECIFICATIONS

Manufactured in:

- Structure and covers manufactured in extruded aluminium, forming removable compartments.
- Cover ends manufactured in polystyrene.

The dimensions are 74 mm depth by 120 mm canal.

Treatments:

- Main body lacquered: **RAL 9010**.
- Covers lacquered in:
 - Red: **RAL 3003**
 - Blue: **RAL 5015**
 - Green: **RAL 6019**
 - Grey: **RAL 7035**
 - White: **RAL 9010**

In compliance with UNE-EN ISO 11197 standard.