



VERTICAL

Vertical bedhead

The **VERTICAL bedhead** is a hospital bedhead designed to respond to the highest aesthetics and functionality requirements, becoming another architectural element of the environment. It is able to integrate the supply of medical gases, electrical, telephone and computer plugs and the call device for nursing assistance, as well as to cast indirect light to the patient's bed.

There are several configurations for these bedheads: they have a standard structure that accommodates one or two Ø38 mm bars for the placement of trays, drawers, monitor or pump supports, etc... Our technical office will adapt the dimensions, colours and final equipment to the specific needs of each project on demand.



TECHNICAL SPECIFICATIONS

Manufactured in:

- Structure and covers in extruded aluminium, forming three independent compartments.
- Translucent covers in polycarbonate.
- Covers ends in ABS.
- Tubes in stainless steel.

Its dimensions are 122 mm depth and 233 mm width.

Finishing:

- Main body lacquered: **RAL 9010**.
- Covers lacquered in the following colours:
 - Red: **RAL 3003**
 - Blue: **RAL 5015**
 - Green: **RAL 6019**
 - Grey: **RAL 7035**
 - White: **RAL 9010**
- The front cover can be customized with different bedhead finishing for different services areas, and the image of the hospital can be incorporated.

In compliance with **UNE-EN ISO 11197** standard.